

# BASIC RADIOGRAPHY X-RAY UNIT (BRXU)



**"CONVERT MOBILE INTO STATIONARY UNIT"**

## DESCRIPTION:


BRXU model is a cost effective, Compact, floor mounted X-ray radiographic unit for all make & types of X-Ray Generators. BRXU is very well suited for nursing homes, diagnostic centre etc with limited amount of space.

## FEATURES:

- Simple sleek, steel fabricated design suitable for small spaces.
- Table top with centre line marking ensures patient centering.
- Simple, smooth & effortless movement of Bucky & column on the entire length of the table with mechanical lock.
- Aluminium railings enables extra Strength with provision of Compression Band.

## SPECIFICATION:

Model		BRXU
Table Top	Material	White sunmica 6mm Thk.
	Table top to Film distance	85 mm
	Table top dimensions	1830 L x 600 W
	Table Top Radiolucent Area	1750 x 540 (mm)
	Load	150 Kg
Bucky Device Unit	Longitudinal Bucky Travel	505 mm (20")
	Locks	Mechanical (Electromagnetic optional)
	Max. Cassette Accommodation	Max 43 cm (17") x 43 cm (17")
	Grid Movement	Reciprocating
	Grid (Standard)	17 1/4 x 18 7/8 60 Lines 6:1 Ratio
	Electrical Requirement	Single Phase 230V +/- 10%. 50/60 Hz
Tube Stand	Tube Column	Integrated Table & Tube Stand
	Vertical Long. Travel	1390 mm (54")
	Weight Balancing Travel	Min 203 mm (8") x 1016 mm (40") max
	Capacity	Upto 40 Kg
Dimension		1830 L x 8800 W x 2030 H (mm)
Weight		Approx Kg.



\*Due to continuous product improvements, specifications & appearance are subject to change without prior notice.

### MANUFACTURER



**SAI  
INDUSTRIES**

SHED NO. D3-12, ANANT INDL. ESTATE PHASE 3,  
BENKAR WASTI, DHAYARI, PUNE—411041.  
MAHARASHTRA (INDIA) (M) +91-9422521699;  
Email: srandive9@gmail.com;  
saiindustries999@gmail.com

### AUTHORISED DISTRIBUTOR