

# X-RAY FIXED HORIZONTAL BUCKY TABLE (FHXT)



## DESCRIPTION:

FHXT model is a cost effective radiographic table to perform fast & simple routine radiography. It is a compact and reliable radiographic table design to meet the requirements of general radiography.

## FEATURES:

- Simple design facilitates quick installation.
- Convenient table top height for easy patient access.
- Table top with centre line marking.
- Simple, smooth & effortless movement of Bucky on the entire length of the table with mechanical lock.
- Sleek, steel fabricated design to take care of the patient weight.
- Aluminium table railings enabling extra strength with provision of Compression Band.



## SPECIFICATION:

Model		FHXT
Table Top	Material	White Sunmica 6mm Thk 324/CG Compact
	Table top to Film distance	85 mm
	Table top dimensions	1830 L x 600 W
	Table Top Radiolucent Area	1750 x 540 (mm)
	Maximum Load	150 Kg
Bucky Device Unit	Reciprocating	Motor & Lever mechanism
	Max. Cassette Accommodation	Max 43 cm (17") x 43 cm (17")
	Electrical requirement	Single Phase 230V +/- 10%. 50/60 Hz
	Grid (Standard)	17 1/4 x 18 7/8 60 Lines 6:1 Ratio
Horizontal Bucky Travel		1050 mm
Locks		Mechanical (Electromagnetic optional)
Dimension		1830 L x 650 W x 655 H (mm)
Weight		Approx 45 Kg.

\*Due to continuous product improvements, specifications & appearance are subject to change without prior notice.

### MANUFACTURER



**SAI  
INDUSTRIES**

SHED NO. D3-12, ANANT INDL. ESTATE PHASE 3,  
BENKAR WASTI, DHAYARI, PUNE—411041.  
MAHARASHTRA (INDIA) (M) +91-9422521699;  
Email: srandive9@gmail.com;  
saiindustries999@gmail.com

### AUTHORISED DISTRIBUTOR