

X-RAY VETERINARY UNIT (VET-X)



DESCRIPTION:

VET-X model is a cost effective radiographic unit for veterinary applications. Designed for small animal veterinarians, VET-X allows maximum imaging coverage with limited installation space.

FEATURES:

- Simple design facilitates quick installation.
- Convenient table top height.
- Table top with centre line marking.
- Simple, smooth & effortless movement of Bucky & column on the entire length of the table with mechanical lock.
- Sleek, steel fabricated design. Aluminium angles for extra strength.


SPECIFICATION:

Model		VET-X
Table Top	Material	White Laminated Sheet 8 mm
	Table top to Film distance	85 mm
	Table top dimensions	1830 L x 600 W
	Table Top Radiolucent Area	1750 x 540 (mm)
	Maximum Load	100 Kg
Bucky Device Unit	Reciprocating	Motor & Lever mechanism
	Max. Cassette Accommodation	Max 43 cm (17") x 43 cm (17")
	Electrical requirement	Single Phase 230V +/- 10%. 50/60 Hz
	Grid (Standard)	17 1/4 x 18 7/8 60 Lines 6:1 Ratio
	Longitudinal Bucky Travel	505 mm (20")
	Locks	Mechanical (Electromagnetic optional)
Tube Stand	Tube Column	Integrated Table & Tube Stand
	Tube Column Height	1856mm (73 ")
	Longitudinal Column Travel	920 mm (36")
	Source to Image Distance	40 " Fixed
Dimension		1830 L x 620 W x 655 H (mm)
Weight		Approx Kg.



*Due to continuous product improvements, specifications & appearance are subject to change without prior notice.

MANUFACTURER



SAI
INDUSTRIES

SHED NO. D3-12, ANANT INDL. ESTATE PHASE 3,
BENKAR WASTI, DHAYARI, PUNE—411041.
MAHARASHTRA (INDIA) (M) +91-9422521699;
Email: srandive9@gmail.com;
saiindustries999@gmail.com

AUTHORISED DISTRIBUTOR