

MOBILE-100

X-RAY GAS SPRING MOBILE STAND (GSXS)



DESCRIPTION:

Mobile Stand GSXS is a very light weight, easy to maneuver with unique gas spring balancing system. It is easily transportable anywhere in the hospital through elevators. It can be folded & transported very safely in hospital wards to take X-Rays in beds without mobilising the patient. It is compatible with 100mA X-ray unit.

FEATURES :

- Gas Spring balancing.
- UP / DN smooth movement.
- Positioning at any FFD dist.
- Safe transportation as Tube Head can be parked to lowest position.
- Wheels & Casters location permits easy positioning for bed / trauma X-Rays
- Foot Lock to immobilize during X-Ray.
- Light weight sleek design.

SPECIFICATION :

- Gas Spring balancing.
- Floor to Focal Spot Dist.
 - Highest – 175 cm
 - Lowest – 35 cm
 - Travel – 140 cm
- Overall Size With Control
 - 1375 (L), 745 (W), 1930 (H).
- Tube Rotation $\pm 180^\circ$
- Tube Swing + 180° / - 45°
- Weight Approx
 - Mobile Stand – 50 kg.

*Due to continuous product improvements, specifications & appearance are subject to change without prior notice.

CONTENTS :

Mobile Stand consists of fabricated

- 1) Base Assembly with wheels, wheel covers, casters & brake assy.
- 2) Vertical Stand with arm, tie bar, gas springs (3 nos.) suitable for balancing 100/60 mA Tube Head & Collimex.



MANUFACTURER



**SAI
INDUSTRIES**

SHED NO. D3-12, ANANT INDL. ESTATE PHASE 3,
BENKAR WASTI, DHAYARI, PUNE—411041.
MAHARASHTRA (INDIA) (M) +91-9422521699;
Email: srandive9@gmail.com;
saiindustries999@gmail.com

AUTHORISED DISTRIBUTOR